



PRODUCT SPECIFICATION



Product: Net Dressing
REF: 1893



1 General description

A large sterile metallized dressing and a latex-free tubular net intended for keeping the dressing in place. A non-adherent layer on both sides of the dressing prevents the dressing from sticking to the wound.

2 Intended use

The dressing is intended to be used as a mechanical barrier and/or for absorption of exudates from large wounds on the head or other parts of the body, e.g. legs.

3 Technical description

Product dimensions

Dressing 20 x 40 cm
Tubular net 7 x 40 cm unstretched

Material

Dressing Absorbent layer of non-woven 100% viscose. Non-adherent layer of polyethylene on both sides of dressing, aluminum metallized on one side.
Tubular net Polyamide and elastomer
Sterile package Paper and transparent PP film.

4 Sterilization

Sterilization method: The metallized dressings are sterilized by irradiation.

5 Composition

Consumer package

<u>Ref#</u>	<u>Product/Material</u>	<u>Pcs</u>
51618171	Metallized dressing	1
51618595	Tubular net	1

Transport package

Consumer packages (pcs):	20
Dimensions (l x w x h) cm:	26 x 20,5 x 17,5
Weight (kg):	0,8

6 Labeling

On product/ consumer package

Cederroth address, telephone- and faxnumber, web-address, E-mail-address
Product name, Reference code
Lot-number, Expiry date, EAN-code, CE

On transport package

Cederroth address, telephone- and faxnumber, web-address
Product name, Reference code
Lot-number, Expiry date, EAN-code, Quantity

7 Storage conditions

The product should be stored in the original package for maximum protection until time of use.
Shelf life: 5 years from production date

Authorization: